

Emergency Trolley



The Task Medical Emergency Trolley is a crucial piece of hospital equipment designed to allow for rapid emergency response and improve organisation and efficiency. The trolley is designed with smooth wheels and brakes for quick manoeuvrability and can easily be transported to different hospital areas, including wards, ICUs, and operating rooms.

Product Features

- Large-capacity with adjustable PC dividers
 - 100mm noise-free castors with brakes and protective bumpers
 - 125 Lightweight aluminium alloy drawers
 - Adjustable system for custom configuration
 - 25kg per drawer, 100kg overall
- Accessories included - Waste bin, wire basket, I.V. Pole and defib shelf

CODE
TKEMT5R



Product Specifications

- Dimensions: L650 x W480 x H1000 mm
- Surface: ABS moulded table with abrasion and corrosion resistance
- Pillars: Aluminium alloy with electrophoresis coating
- Handles: ABS for durability and comfort

Outer Carton



www.constar.com.au

