



CLARITY

LEADING THE WAY IN ULTRASOUND TECHNOLOGY

Our extensive portfolio of ultrasound products is used and trusted by Sonographers throughout Australia. Experience superior quality and performance.

CONS  AR



1300 778 924

Find out more at www.constar.com.au

**Advanced
Formula!**

Clarity Ultrasound Transmission Gel

**EFFORTLESSLY
PUMP AND REFILL
250ML FROM
5L, NO LIFTING
REQUIRED**

Our Clarity Ultrasound Transmission Gel is designed for superior acoustic performance. **Available in both 5L and 250ml bottles, in blue or clear,** it ensures clear, uninterrupted ultrasonography and Doppler applications.

- Suitable for all ultrasound frequencies
- Hypoallergenic, CIT/MIT, and paraben-free
- High viscosity, vacuum-treated
- Water-soluble, oil-free, odourless



**250ML
12 Pack
CLUSG250**



CLUSG5CL

CLUSG5BL

5 Litre

Clarity Gel Warmers

**Maintain Constant
Temperature**

Enhance patient comfort with Clarity Ultrasound Gel Warmers. **Available in single and dual bottle models,** these warmers feature adjustable temperature control and a sleek design, making them an essential addition to any ultrasound setup.

- LED display, adjustable temperature (34-40°C)
- Programmable timer, energy-efficient
- Compact design

**Perfect Transmission
Of Soundwaves**



**Single
CLGW01**



**Double
CLGW02**



1300 778 924

Find out more at www.constar.com.au

Clarity Probe Covers

One Size
Fits Most

Clarity Probe Covers provide reliable protection for your ultrasound probes. **Available in lubricated, non-lubricated, and latex-free options**, these covers ensure hygienic and efficient examinations.

- Single walled
- Ultra smooth surface
- Rolled for convenience



Lubricated
CLPCVL



Non-Lubricated
CLPCV



Latex Free
CLPCVLF



Clarity Freezer Bags

Latex-Free Bands

Our Clarity Sterile Freezer Bags are designed for single-use and come with two latex-free elastic bands. These bags offer a sterile and convenient solution for your medical storage needs.

- 30 x 45cm
- Includes 2 latex-free elastic bands
- Individually packed



30 x 45cm

Individually
Packed



Order Code
CLFBW



1300 778 924

Email us at info@constar.com.au



Enhance Your Ultrasound Experience with These Essential Products

Discover our carefully curated selection of complementary products designed to enhance the efficiency and effectiveness of your ultrasound procedures.

Positioning Pad

A high-density foam wedge designed for optimal patient positioning during female pelvic ultrasound. It features a durable, hospital-grade vinyl cover that is easy to clean and provides a stable, comfortable surface on non-tilting couches.


- Rigid for support
- Soft for comfort
- Easy to clean

150 x 550 x 500mm

Durable & Hygienic

TASK
MEDICAL

Order Code
TKPOSTV



Positioning Pad 45°

Available in a convenient pack of 2, perfect for always having a backup ready when you need it.


Enhanced Stability

255 x 180
x 205mm

258 x 180
x 60mm



Short 45°
TKPOS45




Long 45°
TKPOS45-60

Positioning Pad Kit

A versatile set of 8 pads, each tailored for specific positioning needs, enhancing imaging accuracy.

Precision Positioning



Order Code
TKPOSPAD



1300 778 924

Find out more at www.constar.com.au

Optimize Cardiac Monitoring with Clarity ECG Electrode

Experience accurate and comfortable cardiac assessments with Clarity ECG Electrode. Designed for precision and ease of use in all monitoring applications.


Clarity ECG Electrode

Designed to provide secure and accurate electrical signals for **holter**, **stress test** and **monitoring** purposes while ensuring patient comfort during extended wear. They come conveniently packaged for easy access and storage in healthcare settings.

- **Cover** - One side siliconized PET cover
- **Conductive gel** - Solid gel
- **Carrier** - PE foam 1 mm thick
- **Label** - PE, round Ø24mm, with overprint
- **Sensor** - Ag/AgCl
- **Stud** - Stainless Steel
- **Finger-lift** - A paper sheet for easy electrode removal
- **Latex Free**
- **Shelf life** - 24 months

Convenient
Application

36 x 48mm



Order Code
CLEL3648

Packaged for Quick Access
& Maximum Efficiency



1300 778 924

Email us at info@constar.com.au

