

Basic Dressing Pack



The Task Medical Basic Dressing Packs are essential for maintaining sterile conditions during wound care and medical procedures. Each pack is latex-free and individually sealed to ensure quality and sterility, conforming to EO Sterile standards. Ideal for healthcare settings, these packs enhance patient safety and support sterile protocols.

Contents

- 1 x Sterile Field
- 1 x Dressing Tray
- 1 x Dressing Towel
- 5 x Cotton Balls
- 2 x Blue Forceps
- 1 x Yellow Forceps
- 2 x Gauze Swabs (7.5 x 7.5 cm)



CODE
TKBDP

Product Information

- Individually sealed for sterility
- Packaged in sealed packs of 20
- Two years shelf life
- 160 x 160 mm and 30 mm deep

Packaging

- 20 packs per bag
- 8 bags per carton
- Total pieces per carton - 160

Carton Size

- 50x34x60cm
- G.W 8kg
- N.W 6.5kg

Base Unit



9 342916 020125

Outer Carton



9 342916 020156



www.constar.com.au

