


# ECG Electrode CLEL3648



Clarity ECG electrodes are designed to provide a secure and accurate electrical contact for monitoring purposes while ensuring patient comfort during extended wear. They come conveniently packaged for easy access and storage in healthcare settings.

## Product Information

- **Cover** - one side siliconized PET cover
- **Conductive gel** - solid gel
- **Carrier** - PE foam 1 mm thick
- **Label** - PE, round  $\varnothing 24$  mm, with overprint
- **Sensor** - Ag/AgCl
- **Stud** - Stainless Steel
- **Fingerlift** - A sheet of paper that facilitates removing the electrode from the cover



## Application

- Holter
- Stress Test
- Monitoring

## Packaging

- 5 electrodes/sheet
- 10 sheets/pack
- 40 packs/carton
- 2000/carton

## Carton Size

- 60x20x26cm
- G.W 2.9kg
- N.W 2.5kg

Inner Box



Outer Carton



[www.constar.com.au](http://www.constar.com.au)

