

Freezer Bags



These sterile and transparent freezer bags are designed to provide you with the ultimate assurance when using them.

Product Information

- Single use only
- Transparent sterile packed
- Size 45cm long x 30cm wide
- Latex Free
- Includes two latex free bands



CODE
CLFBW

Packaging

- 100 bags/box
- 20 boxes/carton
- Total pieces/carton - 2000

Carton Size Box Size

- 53x43x51cm ■ 40x10x12cm
- G.W 19.5kg ■ Weight 0.95kg
- N.W 18kg

Inner Box



Outer Carton



www.constar.com.au

