

# Gel Warmer Double



Ensure patient comfort with the Clarity Ultrasound Gel Warmer, available in both single and dual bottle models. Designed with advanced digital temperature control technology, this elegant and compact gel warmer maintains the perfect gel temperature for optimal patient experience. The high-brightness LED display, adjustable temperature control, and programmable timer ensure user-friendly operation and energy efficiency.

## Product Information

- High-brightness LED display
- Adjustable temperature control (34 – 40°C)
- Programmable timer
- Digital temperature control technology
- High energy efficiency
- Elegant and compact design

CODE  
CLGW02



## Packaging

- 1 gel warmer/box

## Box Size

- 29x25x16cm
- Weight 1.20kg



[www.constar.com.au](http://www.constar.com.au)

