

# Latex Free Probe Covers



Clarity Non Latex Probe Covers are a cutting-edge solution designed for endocavity imaging on ultrasound systems. Crafted from non-latex synthetic polyisoprene, this probe cover ensures a versatile fit across various ultrasound transducers, offering healthcare professionals unmatched flexibility during examinations.

## Product Information

- Ultra smooth surface
- Single walled
- Translucent natural colour
- Length : 230 mm
- Width : 53 mm ± 2 mm
- Diameter : 33 mm ± 1.5 mm
- Thickness : 0.085 ± 0.5 mm
- Latex free



## Packaging

- 100 probe covers/box
- 20 boxes/carton
- 2000 pieces/carton

## Carton Size

- 49x39x32cm
- G.W 8kg
- N.W 8.64kg

## Box Size

- 25.6x13x7cm
- Weight 0.40kg

Inner Box



Outer Carton



[www.constar.com.au](http://www.constar.com.au)

