


Fenestrated Drape 50cm x 68cm



The Task Medical Fenestrated Drape is ideal for various surgical and therapeutic procedures, helping protect patients and medical staff from bacterial transfer.

Product Information

- Peel pack
- Antistatic
- 4 ply poly-lined
- Individually sterile packed
- Folded design that is easy to open
- 10cm diameter fenestration
- Non-woven fabric
- Single use only
- EO sterilised



Packaging

- 20 pouches/box
- 18 boxes/carton
- Total pieces/carton - 360

Carton Size

- 54x32x45cm
- G.W 11kg
- N.W 9kg

Box Size

- 20x15x17cm
- Weight 0.50kg

Base Unit



Outer Carton



www.constar.com.au

