

# Isolation Gown



Task Medical Disposable SMS Isolation Gowns provide an effective hygienic barrier against contamination, protecting the body and clothing from splashes, soiling, dust, and dirt. Made from superior 25gsm SMS non-woven fabric, these gowns are lightweight yet durable, ensuring comfort and safety for the wearer. Ideal for use in medical settings, they feature knitted cuffs, waist and neck ties.

## Product Information

- Made from 25gsm SMS non-woven fabric (two Spunbond layers and one Meltblown layer)
- Light blue colour
- Knitted cuffs for a secure fit
- Disposable
- Neck and waist ties for easy adjustment
- Provides protection for arms, torso, and upper to mid legs
- 140 × 120 cm, regular fit



CODE  
TKGOWNISOB

## Packaging

- 10 sheets per pack
- 5 packs per carton
- Total pieces per carton - 50

## Carton Size

- 45x27x32cm
- G.W 5.7kg
- N.W 5.2kg

Outer Carton



[www.constar.com.au](http://www.constar.com.au)

