

# Isolation Gown



Task Medical Disposable SMS Isolation Gowns provide an effective hygienic barrier against contamination, protecting the body and clothing from splashes, soiling, dust, and dirt. Made from a superior 25gsm SMS non-woven fabric, these gowns are lightweight yet durable, ensuring comfort and safety for the wearer. Ideal for use in medical settings, they feature knitted cuffs, waist and neck ties.

## Product Information

- Made from 25gsm SMS non-woven fabric (two Spunbond layers and one Meltblown layer)
- Light blue color
- Knitted cuffs for a secure fit
- Disposable
- Neck and waist ties for easy adjustment
- Provides protection for arms, torso, and upper to mid legs
- 140 × 120 cm, regular fit



CODE  
TKGOWNISOB

## Packaging

- 10 sheets per pack
- 5 packs per carton
- Total pieces per carton - 50

## Carton Size

- 45x27x32cm
- G.W 5.7kg
- N.W 5.2kg

Outer Carton



[www.constar.com.au](http://www.constar.com.au)

