

Mixing Cannula 14cm



The Task Medical Mixing Cannula is a semi-rigid, sterile tube designed to facilitate accurate, contamination-free mixing of liquids or solutions in healthcare, laboratory, and research environments.

Product Information


- Peel pack
- Precise reach and manoeuvrability
- Compatible with various solutions or fluids
- User-friendly design with a tapered tip
- Secure attachment mechanism
- Ideal for sterile mixing and chemical applications
- Complies with medical industry standards
- Single-use for infection control
- Designed to be attached to standard luer slip syringes
- EO sterilised

Packaging

- 100 pieces/bag

Carton Size

- 41x28x18.5cm
- G.W 3.05kg
- N.W 2.55kg



CODE
TKMIX

14 cm long

Base Unit



Outer Carton



www.constar.com.au

