

# Premium Disposable Fitted Bed Sheet



Task Medical Premium Fitted Bed Sheets snugly fit most medical examination couches, offering a protective barrier between patients and furniture. Made from 40gsm SMS non-woven fabric, these sheets ensure hygiene and cleanliness in medical settings.

## Product Information

- 40gsm SMS non-woven fabric
- Fitted corners, elastic head and foot
- Light blue, soft, and comfortable
- Disposable, fluid resistant, tear-resistant
- 200 x 70 x 10 cm

## Packaging

- 10 sheets per pack
- 10 packs per carton
- Total pieces per carton - 100

## Carton Size

- 50x30x50cm
- G.W 8.7kg
- N.W 8kg



Inner Box



Outer Carton



[www.constar.com.au](http://www.constar.com.au)

