

# Premium Disposable Modesty Sheets



The Task Medical Premium Modesty Sheet is made from high-quality 40gsm SMS non-woven fabric, ensuring exceptional protection against contamination. Designed for patient privacy and hygiene, it is perfect for covering patients or as an undersheet. Its smaller size makes it versatile for various medical applications.

## Product Information

- 40gsm SMS non-woven fabric
- Light blue, soft, and comfortable
- Disposable, fluid resistant, tear-resistant
- Prevents liquid leakage
- Higher elasticity
- Breathable
- Excellent barrier
- Protects patients and furniture
- 100 x 70 cm



## Packaging

- 10 sheets per pack
- 10 packs per carton
- Total pieces per carton - 100

## Carton Size

- 45x27x20cm
- G.W 6kg
- N.W 5kg

Outer Carton



[www.constar.com.au](http://www.constar.com.au)

