


# Premium Disposable Pillow Cases



Task Medical Premium Disposable Pillow Cases are meticulously designed to cover most medical pillows, offering a protective barrier for patients and furniture surfaces. Made from superior 40gsm SMS non-woven fabric, these light blue pillow cases provide exceptional hygiene and cleanliness. Ideal for maintaining a pristine environment, they shield examination tables from contamination caused by body oils, perspiration, fluids, and cosmetics.

## Product Information

- 40gsm SMS non-woven fabric
- Sealed on three sides, open end with tuck flap
- Light blue color
- Soft and comfortable
- Disposable, fluid resistant
- Prevents liquid leakage
- High elasticity, tear-resistant
- Breathable and effective barrier
- 76 x 50 cm



CODE  
TKPRPLCS-LB

## Packaging

- 10 sheets per pack
- 20 packs per carton
- Total pieces per carton - 200

## Carton Size

- 55x27x28cm
- G.W 6.5kg
- N.W 5.5kg

Inner Box



Outer Carton



[www.constar.com.au](http://www.constar.com.au)

