

SMS Isolation Gown



Task Medical Disposable SMS Isolation Gowns provide an effective hygienic barrier against contamination, protecting the body and clothing from splashes, soiling, dust, and dirt. Made from superior 25gsm SMS non-woven fabric, these gowns ensure comfort and safety with their lightweight design. Ideal for medical settings, they feature knitted cuffs, waist and neck ties.

Product Information

- Made from 25gsm SMS non-woven fabric (two Spunbond layers and one Meltblown layer)
- Light blue colour
- Knitted cuffs for a secure fit
- Disposable
- Neck and waist ties for easy adjustment
- Provides protection for arms, torso, and upper to mid legs
- 140 × 120 cm, regular fit



Packaging

- 10 sheets per pack
- 5 packs per carton
- Total pieces per carton - 50

Carton Size

- 29x27x32cm
- G.W 3.4kg
- N.W 3kg

Outer Carton



www.constar.com.au

