

# Shoe Cover



The Task Non-Skid Non-Woven Disposable Shoe Cover in blue offers reliable protection in various medical and industrial settings. Featuring anti-skid painting on the sole and ultrasonic seams, these shoe covers ensure safety and durability.

## Product Information

- Non-skid sole for secure footing
- Impervious material suitable for heavy fluid applications
- Latex-free and hypoallergenic
- Recyclable polypropylene construction
- 16 inches, one size fits all

## Packaging

- 100 covers per pack
- 10 packs per carton
- Total pieces per carton - 1000



## Carton Size

- 33x25x28cm
- G.W 6kg
- N.W 5.6kg

Base Unit



Outer Carton



[www.constar.com.au](http://www.constar.com.au)

