

Tourniquet Multicolor

TASK MEDICAL

The Task Medical Multicoloured Tourniquet is designed for daily use in hospitals, medical practices, and by paramedics. Crafted meticulously for reliability and comfort, this tourniquet ensures ease of use and patient safety.

Product Information

- Skin-friendly, flexible, soft ribbon
- 25 mm wide
- Sturdy safety snap closure
- One-hand operation for convenience
- Replaceable and washable ribbon
- Available with or without quick-release lock
- Metal and latex-free
- No squeezing or pinching of the skin
- Made in Germany



Inner Box



Outer Carton



www.constar.com.au

