

# Tube Connector



Task Medical Sterile Tube Connector helps to facilitate secure and hygienic connections between tubing in clinical and laboratory settings.

## Product Information

- Peel pack
- Size medium 7 - 10mm
- EO sterilised
- Versatile for different tubing sizes
- Secure, leak-proof connection between tubes
- Easy handling and installation
- Transparency for visual inspection
- Complies with medical industry standards
- Single-use for infection control
- Ideal for various medical applications

CODE  
TKTCM



## Packaging

- 50 pieces/bag
- 2 bags/carton
- 100 pieces/carton

## Carton Size

- 34x16x15cm
- G.W 0.55kg
- N.W 0.35kg

Base Unit



Outer Carton



[www.constar.com.au](http://www.constar.com.au)

