

# Latex Examination Gloves - Extra Large

# assurity latex

Assurity latex examination gloves in extra large are designed for comfort and flexibility. Made from natural rubber latex, these powder-free gloves provide a roomy fit, ideal for long wear in healthcare environments.

## Product Information

- Premium quality latex
- Powder-free and micro-textured
- Comfortable and ambidextrous
- Beaded cuff for easy donning
- Non-sterile, single use
- AQL 1.5
- Eco-packed (100 gloves per box)



## Packaging

- 100 gloves per box
- 10 boxes/carton
- 1000 pieces/carton

## Carton Size

- 36.5x26x24.4cm
- G.W 7.7kg
- N.W 6.7kg

## Box Size

- 23.5x12.5x70cm
- Weight 0.70kg

Base Unit



Outer Carton



[www.constar.com.au](http://www.constar.com.au)



# Latex Examination Gloves

# assurity latex

The Assurity latex examination gloves range provides healthcare professionals with quality and protection. Available in five sizes, these powder-free gloves offer excellent flexibility and grip. Each box is colour-coded for easy identification and tested to high international standards.

## Product Information

- Box of 100 gloves
- Made from natural rubber latex
- Powder-free, ambidextrous, micro-textured
- Palm thickness 4.5mil, finger thickness 6.0mil
- Beaded cuff, non-sterile, single use
- ASTM D3578 | ARTG ID 392396
- Conform to EN374-1, EN374-5, EN455-1,2,3,4, EN21420 standards
- Certified ISO9001:2015, ISO13485:2016



Extra Small  
ASLEGNXS



Small  
ASLNEGNS



Medium  
ASLEGNM



Large  
ASLEGNL



Extra Large  
ASLEGNXL



[www.constar.com.au](http://www.constar.com.au)

