

Latex Examination Gloves - Extra Small

assurity latex

Assurity latex examination gloves in extra small size are designed for healthcare professionals, offering reliable barrier protection. Made from natural rubber latex, these powder-free gloves ensure comfort and dexterity.

Product Information

- Premium quality latex
- Powder-free and micro-textured
- Comfortable and ambidextrous
- Beaded cuff for easy donning
- Non-sterile, single use
- AQL 1.5
- Eco-packed (100 gloves per box)



Packaging

- 100 gloves per box
- 10 boxes/carton
- 1000 pieces/carton

Carton Size

- 36.5x26x24.4cm
- G.W 5.9kg
- N.W 4.9kg

Box Size

- 23.5x12.5x70cm
- Weight 0.70kg

Base Unit



Outer Carton



www.constar.com.au



Latex Examination Gloves

assurity latex

The Assurity latex examination gloves range provides healthcare professionals with quality and protection. Available in five sizes, these powder-free gloves offer excellent flexibility and grip. Each box is colour-coded for easy identification and tested to high international standards.

Product Information

- Box of 100 gloves
- Made from natural rubber latex
- Powder-free, ambidextrous, micro-textured
- Palm thickness 4.5mil, finger thickness 6.0mil
- Beaded cuff, non-sterile, single use
- ASTM D3578 | ARTG ID 392396
- Conform to EN374-1, EN374-5, EN455-1,2,3,4, EN21420 standards
- Certified ISO9001:2015, ISO13485:2016



Extra Small
ASLEGNXS



Small
ASLNEGNS



Medium
ASLEGNM



Large
ASLEGNL



Extra Large
ASLEGNXL



www.constar.com.au

