

Latex Examination Gloves - Medium

assurity latex

Assurity latex examination gloves in medium offer a balanced combination of protection and dexterity. Made from natural rubber latex, these powder-free gloves provide a secure fit for healthcare professionals.

Product Information

- Premium quality latex
- Powder-free and micro-textured
- Comfortable and ambidextrous
- Beaded cuff for easy donning
- Non-sterile, single use
- AQL 1.5
- Eco-packed (100 gloves per box)



Packaging

- 100 gloves per box
- 10 boxes/carton
- 1000 pieces/carton

Carton Size

- 36.5x26x24.4cm
- G.W 6.7kg
- N.W 5.7kg

Box Size

- 23.5x12.5x70cm
- Weight 0.70kg

Base Unit



Outer Carton



www.constar.com.au



Latex Examination Gloves

assurity latex

The Assurity latex examination gloves range provides healthcare professionals with quality and protection. Available in five sizes, these powder-free gloves offer excellent flexibility and grip. Each box is colour-coded for easy identification and tested to high international standards.

Product Information

- Box of 100 gloves
- Made from natural rubber latex
- Powder-free, ambidextrous, micro-textured
- Palm thickness 4.5mil, finger thickness 6.0mil
- Beaded cuff, non-sterile, single use
- ASTM D3578 | ARTG ID 392396
- Conform to EN374-1, EN374-5, EN455-1,2,3,4, EN21420 standards
- Certified ISO9001:2015, ISO13485:2016



Extra Small
ASLEGNXS



Small
ASLNEGNS



Medium
ASLEGNM



Large
ASLEGNL



Extra Large
ASLEGNXL



www.constar.com.au

