

Nitrile Examination Gloves - Extra Large

assurity nitrile

Assurity nitrile examination gloves in extra large are designed for maximum comfort and protection. Constructed from 100% nitrile, these powder-free gloves offer a spacious fit, ideal for prolonged wear in demanding environments.

Product Information

- Premium quality nitrile
- Powder-free and latex-free
- Micro-textured for improved grip
- Ambidextrous and durable
- Chemotherapy tested
- AQL 1.5
- Eco-packed (200 gloves per box)



Packaging

- 200 gloves per box
- 10 boxes/carton
- 2000 pieces/carton

Carton Size

- 43x25.5x24.8cm
- G.W 9.2kg
- N.W 8.2kg

Box Size

- 24x12.5x85cm
- Weight 0.80kg

Base Unit



Outer Carton



www.constar.com.au



Nitrile Examination Gloves

assurity nitrile

The Assurity nitrile examination gloves range provides healthcare professionals with reliable protection and comfort. Available in five sizes, these gloves mimic the feel of latex while reducing hand fatigue. Their micro-textured finish enhances grip, making them suitable for a variety of tasks. Each box is colour-coded for easy identification and rigorously tested to meet international standards.

Product Information

- Box of 200 gloves
- 100% nitrile material
- Powder-free, latex-free
- Chemotherapy tested, micro-textured finish
- 3.2mil nitrile / Finger thickness 4.0mil, non-sterile, single use
- Beaded cuff for easy donning
- ASTM D6319 & ASTM D6978 | ARTG ID 392394
- Conform to EN374-1, EN374-5, EN455-1,2,3,4, EN21420 standards
- Certified ISO9001:2015, ISO13485:2016



Extra Small
ASNEGBXS



Small
ASNEGBS



Medium
ASNEGBM



Large
ASNEGBL



Extra Large
ASNEGBXL



www.constar.com.au

