

Nitrile Examination Gloves - Medium

assurity
nitrile

Assurity nitrile examination gloves in medium deliver a perfect balance of protection and dexterity. Constructed from 100% nitrile, these powder-free gloves ensure a comfortable fit for various healthcare tasks.

Product Information

- Premium quality nitrile
- Powder-free and latex-free
- Micro-textured for improved grip
- Ambidextrous and durable
- Chemotherapy tested
- AQL 1.5
- Eco-packed (200 gloves per box)



Packaging

- 200 gloves per box
- 10 boxes/carton
- 2000 pieces/carton

Carton Size

- 43x25.5x24.8cm
- G.W 8.0kg
- N.W 7.0kg

Box Size

- 24x12.5x85cm
- Weight 0.80kg

Inner Box



Outer Carton



www.constar.com.au



Nitrile Examination Gloves

assurity nitrile

The Assurity nitrile examination gloves range offers healthcare professionals reliable protection and comfort. Available in five sizes, these gloves mimic the feel of latex while reducing hand fatigue. The micro-textured finish provides enhanced grip, making them suitable for various tasks. Each box is colour-coded for easy identification and tested to meet international standards.

Product Information

- Box of 200 gloves
- 100% nitrile material
- Powder-free, latex-free
- Chemotherapy tested, micro-textured finish
- 3.2mil nitrile / Finger thickness 4.0mil, non-sterile, single use
- Beaded cuff for easy donning
- ASTM D6319 & ASTM D6978 | ARTG ID 392394
- Conform to EN374-1, EN374-5, EN455-1,2,3,4, EN21420 standards
- Certified ISO9001:2015, ISO13485:2016



Extra Small
ASNEGBXS



Small
ASNEGBS



Medium
ASNEGBM



Large
ASNEGBL



Extra Large
ASNEGBXL



www.constar.com.au

