

# Nitrile Examination Gloves - Small

# assurity nitrile

Assurity nitrile examination gloves in small offer reliable barrier protection and enhanced sensitivity. Crafted from 100% nitrile, these powder-free gloves provide a secure fit for users requiring agility and comfort.

## Product Information

- Premium quality nitrile
- Powder-free and latex-free
- Micro-textured for improved grip
- Ambidextrous and durable
- Chemotherapy tested
- AQL 1.5
- Eco-packed (200 gloves per box)



## Packaging

- 200 gloves per box
- 10 boxes/carton
- 2000 pieces/carton

## Carton Size

- 43x25.5x24.8cm
- G.W 7.6kg
- N.W 6.6kg

## Box Size

- 24x12.5x85cm
- Weight 0.80kg

Base Unit



Outer Carton



[www.constar.com.au](http://www.constar.com.au)



# Nitrile Examination Gloves

# assurity nitrile

The Assurity nitrile examination gloves range provides healthcare professionals with reliable protection and comfort. Available in five sizes, these gloves mimic the feel of latex while reducing hand fatigue. Their micro-textured finish enhances grip, making them suitable for a variety of tasks. Each box is colour-coded for easy identification and rigorously tested to meet international standards.

## Product Information

- Box of 200 gloves
- 100% nitrile material
- Powder-free, latex-free
- Chemotherapy tested, micro-textured finish
- 3.2mil nitrile / Finger thickness 4.0mil, non-sterile, single use
- Beaded cuff for easy donning
- ASTM D6319 & ASTM D6978 | ARTG ID 392394
- Conform to EN374-1, EN374-5, EN455-1,2,3,4, EN21420 standards
- Certified ISO9001:2015, ISO13485:2016



Extra Small  
ASNEGBXS



Small  
ASNEGBS



Medium  
ASNEGBM



Large  
ASNEGBL



Extra Large  
ASNEGBXL



[www.constar.com.au](http://www.constar.com.au)

