

EasyRoller Wheelchair

MRI SAFE

TASK

MEDICAL

The EasyRoller MRI Safe Wheelchair is made from 100% recyclable thermoplastics, designed for safe and comfortable patient transport in MRI environments. This wheelchair has no metal parts, ensuring safety and compatibility with MRI systems.

Product Features

- No metal parts
- Recyclable plastic
- Ergonomic design
- Easy manoeuvring
- Maintenance-free

Specifications

- Chair dimensions 67 cm x 82.5 cm x 105 cm
- Weight 20kg
- Seat Dimensions 50 cm x 42 cm x 52 cm
- SWL: 150kg



Carton



www.constar.com.au

