

Fixed Height Stretcher Trolley

MRI SAFE

TASK

MEDICAL

This non-magnetic stainless steel stretcher is tested with the SIEMENS 3.0 Tesla system and is ideal for NMR, CT rooms, and veterinary hospitals. It meets stringent standards with CE and MR compatibility certifications.

**Compatible with 1.5T and 3.0T MRI systems, including
SIEMENS, GE, PHILIPS, and more.**

Product Features

- Stainless steel construction
- Drop-down side rails
- 5" lockable castors (2 with brakes, 2 without)
- Adjustable headrest
- 2" pad included

Product Specifications

- Length: 76.77" (195cm)
- Width: 23.6" (60cm)
- Height: 31.5" (80cm)
- Max Load: 159kgs
- Back Adjustment: 0-60°



Carton



www.constar.com.au

