

Hydraulic Stretcher Trolley - Flat Top

MRI SAFE

TASK MEDICAL

The Task Medical MRI Conditional Adjustable Height Stretcher is a non-magnetic, height-adjustable stretcher, fully tested to 3T, ideal for use in MRI rooms. Its robust design ensures patient safety and ease of use in MRI environments, with smooth adjustability for various patient needs.

Product Features

- Non-magnetic materials
- Height adjustable from 45cm to 85cm
- Fully tested to 3T
- Easy to clean
- Durable for long-term use

Specifications

- Size: L195 x W69 x H45-85CM
- Max Capacity: 200kgs
- Max Height: 85cm
- Min Height: 45cm



Carton



www.constar.com.au

