

# Hydraulic Stretcher Trolley

MRI SAFE

**TASK**  
MEDICAL

The Task Medical MRI Safe Hydraulic Height Adjustable Stretcher prioritises patient safety and comfort, with a fully tested design that ensures reliable performance in MRI environments. With height adjustability and easy manoeuvrability, this stretcher is an essential tool in MRI and clinical care.

## Product Features

- Aluminium alloy structure
- Adjustable backrest to 70°
- IV pole for added convenience
- Hydraulic dual pedal pump for height adjustment
- Trendelenburg supplied as standard
- Non-magnetic castors for MRI safety

## Product Specifications

- Length: 190 cm
- Width: 69 cm
- Height Adjustable: 65 - 102 cm
- Backrest: 0 -70° ( $\pm 2^\circ$ )
- Trendelenburg: 0 -12° ( $\pm 2^\circ$ )
- Max Weight: 200kg



Carton



[www.constar.com.au](http://www.constar.com.au)

