

# Patient Transport Board

MRI SAFE

# TASK

MEDICAL

Facilitate smooth and safe patient transfers with the Task Medical MRI/MR Safe Patient Transport Board. Made from high-density polyethylene, it ensures ease of use in MRI environments with six carrying handles for effective and safe transfers.

## Product Features

- High-density polyethylene
- MRI safe and latex-free
- Six carrying handles
- Easy to clean
- Safe for various patients

## Specifications

- Dimensions: 120cm x 58cm x 0.5cm
- Weight Capacity: 250kgs



Carton



[www.constar.com.au](http://www.constar.com.au)

