

3-Function Hospital Electric Bed

TASK
MEDICAL

The Task Medical Three-Function Electric Bed provides comfort and safety with a remote-controlled adjustable backrest, leg rest, and height. Durable side rails, corner bumpers, and removable head and footboards ensure patient security and comfort.

Product Features

- Electric and Manual CPR for quick adjustments
- Four quiet, durable motors with built-in control buttons
- Side rails with slow-down design, corner bumpers, and removable head/footboard
- Liquid angle indicators, bedside card display, infusion holes, and drainage hooks
- Wood Grain finish
- Nurse control/remote control



Product Specifications

- Dimensions: L2250 x W1040 x H450-750 mm
- Bed Platform: 1980 x 900 mm
- Backrest: 0° - $70^{\circ}\pm 5^{\circ}$
- Leg rest: 0° - $40^{\circ}\pm 5^{\circ}$
- Height: 490 - 790 ± 10 mm
- Trendelenburg: 0° - $12^{\circ}\pm 3^{\circ}$
- Brakes: Central brake

Outer Carton



www.constar.com.au

