

5-Function Hospital Electric Bed (Build-in-Button)

TASK

MEDICAL

The Task Medical Five-Function Electric Bed provides ultimate comfort and safety for patients in hospital or clinical settings. It features fully adjustable backrest, leg rest, height, and Trendelenburg/Reverse Trendelenburg positioning. Ergonomic side rails, a removable headboard and footboard, a liquid angle indicator, electric CPR functionality, corner bumpers, and four motors ensure safety, convenience, and easy operation.

Product Features

- Electric and Manual CPR for quick adjustments
- Four quiet, durable motors with built-in control buttons
- Side rails with slow-down design, corner bumpers, and removable head/footboard
- Liquid angle indicators, bedside card display, infusion holes, and drainage hooks
- Wood Grain finish
- Nurse control/remote control



Product Specifications

- Dimensions: L2250 x W1040 x H490-790 mm
- Bed Platform: 1980 x 900 mm
- Backrest: 0°-70°±5°
- Leg rest: 0°-40°±5°
- Height: 490-790±10 mm
- Trendelenburg: 0°-12°±3°
- Brakes: Central brake

Outer Carton



www.constar.com.au

