

## 5-Function Hospital Electric Bed

# TASK MEDICAL

The Task Medical 5-Function Hospital Electric Bed offers maximum comfort and safety with fully adjustable settings for backrest, leg rest, and height for patients and caregivers. Features like slow-down side rails, corner bumpers, and a quiet motor to ensure patient security and a peaceful environment. The bed also includes removable head and footboards and optional extensions for flexible patient care.

### Product Features

- Manual and Electric CPR for emergency adjustments
- Four high-quality motors, controlled via remote
- Damping side rails, removable head/footboard, and corner bumpers
- Blue angle indicators, infusion holes, and drainage hooks
- Central brake, 125mm castors, remote control
- Wood Grain finish
- Auto-contour and calf raising options
- Manual CPR function
- Adjustable I.V Pole



### Product Specifications

- Dimensions: L2250 x W1040 x H490-790 mm
- Bed Platform: 1980 x 900 mm
- Backrest:  $0^{\circ}-70^{\circ}\pm 5^{\circ}$
- Leg rest:  $0^{\circ}-40^{\circ}\pm 5^{\circ}$
- Height:  $490-790\pm 10$  mm
- Trendelenburg & Reverse Trendelenburg:  $0^{\circ}-12^{\circ}\pm 3^{\circ}$

Outer Carton



[www.constar.com.au](http://www.constar.com.au)

