

# Anaesthesia Trolley

# TASK MEDICAL

The Task Medical Anaesthesia Trolley allows anaesthetists to focus on patient care by having a trolley that can be easily transported between operating rooms, ICUs, and emergency areas. Designed with easy-to-clean, antimicrobial surfaces to minimize contamination and ensures quick access to critical medications and instruments during procedures.

## Product Features

- Spacious and organised for anaesthesia supplies
- Four  $\phi 125$ mm castors, two with brakes
- Loading capacity, 25kg per drawer, 100kg total
- Adjustable shelf customisable for various storage needs
- PC Divider, durable, chemical-resistant, and heat/cold resistant
- Waste bucket and basket wire included

## Product Specifications

- Dimensions: L650 x W480 x H993 mm
- Surface: ABS moulded table, abrasion and corrosion-resistant
- Pillars: Aluminium alloy with electrophoresis coating
- Board: Aluminium composite for strength
- Drawers: Zinc-coated steel with detachable PC dividers
- Handles: Durable ABS for easy grip

CODE  
TKAT5BLU



Outer Carton



[www.constar.com.au](http://www.constar.com.au)

