

# Cabinet Trolley



The Task Medical Cabinet Trolley is a durable, mobile storage solution designed for medical environments. Made with cold-rolled steel, PP, ABS, and aluminium alloy, it offers secure and organised storage. Features include a sliding door, adjustable medicine basket, and a robust handle for easy movement.

CODE  
TKCABT

## Product Features

- Sliding Door, PP material, easy to clean and disassemble
- Medicine Basket, transparent PC dividers, adjustable for easy organisation
- Secure cam lock for safety
- Four double castors, two with brakes for stability
- Ergonomic aluminium alloy handle integrated for comfort

## Product Specifications

- Dimensions: L822 x W658 x H1635 mm
- Cabinet Body: Cold-rolled steel with electrostatic spray for durability
- Sliding Door: PP material with ABS sliding grooves
- Castors: Four double castors, two with brakes
- Handle: Aluminium alloy integrated molding



Outer Carton



[www.constar.com.au](http://www.constar.com.au)

