

Infusion Trolley

TASK MEDICAL

The Task Medical Infusion Trolley keeps all necessary IV equipment in one place with large classified storage drawers. Designed to allow healthcare professionals to have quick access to IV supplies and medications to respond to patient needs.

Product Features

- Classified with extra space underneath
- Four ϕ 125mm castors, noise-free with two brakes
- Loading capacity, 25kg per drawer, 100kg total
- Innovation divider, durable polycarbonate for insulation
- Integrated accessories system, adjustable for easy reconfiguration

Product Specifications

- Dimensions: L650 x W480 x H888 mm
- Surface: ABS moulded with plastic plate, abrasion and corrosion-resistant
- Pillar: Aluminium alloy with electrophoresis coating
- Board: Aluminium composite for strength
- Drawers: Zinc-coated steel with PC dividers, detachable for easy cleaning
- Handle: Durable ABS for firm grip

CODE
TKITBLU



Outer Carton



www.constar.com.au

