

Medicine Trolley

TASK MEDICAL

The Task Medical Medicine Trolley is an essential piece of hospital and healthcare equipment used to store, organize, and transport medications safely. This trolley can easily be transported across hospital wards, clinics, and care homes whilst minimising medication errors and allows for immediate access to essential medications. Designed with infection control in mind, the surfaces are easy to clean to prevent cross-contamination.

Product Features

- 8 tray layers, holds 48 medicine boxes for efficient dispensing
- Sliding covers, easy monitoring of medicines
- Prevents spills with secure pull-out trays
- Adjustable partitions, categorise by time (e.g., morning, evening)
- Colour-coding: Blue for undispensed, red for dispensed (medicines/beds)
- 4 noise-free $\phi 125$ mm castors with brakes
Loading capacity 25kg per drawer, 100kg total

CODE
TKMEDTBLU



Product Specifications

- Dimensions: L650 x W480 x H965 mm
- Surface: ABS with abrasion and corrosion resistance, plus a rim to secure items
- Pillar: Aluminium alloy with an electrophoresis coating for durability
- Medicine Boxes: Detachable trays and boxes for hygienic and long-lasting use

Outer Carton



www.constar.com.au

