

Premium Emergency Trolley

TASK MEDICAL

The Task Medical Premium Emergency Trolley is crafted to enhance efficiency, safety, and mobility in critical situations for healthcare professionals.

Product Features

- Ergonomic, cloud-shaped ABS base with a seamless, easy-to-clean design
- Pop-up side table for additional workspace
- Five aluminium drawers: 2x 110mm, 2x 150mm, 1x 210mm, featuring soft-close rails, adjustable dividers, and waterproof label slots
- Embedded seal lock with recordable numbers for secure access
- Four 125mm castors with brakes for stability and smooth movement
- Includes oxygen cylinder and power storage slots for convenience
- Waste bin, I.V. Pole, wire basket and defib shelf included



CODE
TKEMTAL5GRY

Product Specifications

- Dimensions: L700 x W532 x H1007 mm
- Material: High-quality ABS, aluminium alloy columns with oxidation treatment
- Castors: 125mm single-sided with brakes
- Handles: Ergonomic ABS with non-slip grip
- Drawers: Adjustable PC dividers, soft-close rails

Outer Carton



www.constar.com.au

