

Treatment Trolley 1 Drawer

TASK MEDICAL

The Task Medical Hospital Treatment Trolley keeps medical tools, dressings, syringes, and medications neatly arranged to ensure quick access to essential treatment supplies. Made with high-strength ABS, PC, and aluminium alloy, it offers practical storage, secure locking systems, and smooth mobility.

Product Features

- PC dividers for safe, organised, and chemical-resistant storage
- Adjustable and removable accessory slots for flexible organisation
- Large drawers with a lower partition for additional storage
- Four noise-free $\phi 125$ mm castors, with two brakes for stability
- Drawers support up to 25kg; total capacity is 100kg
- Waste bin included

CODE
TKIDTRTBLU



Product Specifications

- Dimensions: L650 x W480 x H888 mm
- Surface: ABS with abrasion and corrosion resistance, plus a rim to secure items
- Pillar: Aluminium alloy with an electrophoresis coating for durability
- Drawers: Zinc-coated steel with PC dividers for easy cleaning and replacement
- Handles: ABS material for long-lasting performance

Outer Carton



www.constar.com.au

