

## Hi-Lo Exam Couch 2 Sections 1 Motor - Navy Blue

# TASK

MEDICAL

The Task Medical Hi-Lo Examination Couch is a 2 Section, motorised examination bed designed for efficiency and patient comfort. With a 70 cm wide surface, it accommodates various patient sizes, making it a reliable choice for healthcare professionals.

### Product Features

- Motorised height adjustment for seamless positioning
- 2-section design with adjustable backrest (-33° to 90°)
- Generous 70 cm width for enhanced patient comfort
- Robust frame for long-term durability
- Lockable wheels for secure manoeuvrability
- Seamless, hygienic navy blue upholstery
- Detachable headrest for versatile procedures

CODE  
TK2SECNB



### Product Specifications

- Flat size: L 1980 mm x W 700 mm
- Height adjustability: 480 – 900 mm (electric actuator with foot switch)
- Safe working load: 250 kg

Outer Carton



[www.constar.com.au](http://www.constar.com.au)

