

## Hi-Lo Exam Couch 3 Sections 1 Motor - Black Edition

# TASK

MEDICAL

The Task Medical Hi-Lo Examination Couch is a 3 Section, motorised examination bed designed for maximum versatility and patient comfort. With a 70 cm wide surface, it accommodates various patient sizes and medical procedures, making it an excellent choice for healthcare professionals. Now available in Black Edition, it offers a sleek and professional aesthetic for modern medical settings.

### Product Features

- Motorised height adjustment for seamless positioning
- 3-section design with adjustable backrest (-33° to 90°) and foot section (-24° to 41°)
- Generous 70 cm width for enhanced patient comfort
- Robust frame for long-term durability
- Retractable castors and levelling feet for stability and manoeuvrability
- Face hole for additional procedure adaptability
- Seamless, hygienic black upholstery for a sleek and professional look

CODE  
TK3SECIMBLKE



### Product Specifications

- Flat size: L 1980 mm x W 700 mm
- Height adjustability: 480 – 900 mm (electric actuator with foot switch)
- Safe working load: 250 kg

Outer Carton



[www.constar.com.au](http://www.constar.com.au)

