

**Hi-Lo Exam Couch
3 Sections 1 Motor
- Black**

TASK
MEDICAL

The Task Medical Hi-Lo Examination Couch is a 3 Section, motorised examination bed designed for maximum versatility and patient comfort. With a 70 cm wide surface, it accommodates various patient sizes and medical procedures, making it an excellent choice for healthcare professionals.

Product Features

- Motorised height adjustment for seamless positioning
- 3-section design with adjustable backrest (-33° to 90°) and foot section (-24° to 41°)
- Generous 70 cm width for enhanced patient comfort
- Robust frame for long-term durability
- Retractable castors and levelling feet for stability and manoeuvrability
- Face hole for additional procedure adaptability
- Seamless, hygienic upholstery also available in Grey, and Black

CODE
TK3SECIMBLK



Product Specifications

- Flat size: L 1980 mm x W 700 mm
- Height adjustability: 480 – 900 mm (electric actuator with foot switch)
- Safe working load: 250 kg

Outer Carton



www.constar.com.au

