

## Hi-Lo Exam Couch 3 Sections 1 Motor Foot Control Bar - Navy Blue

# TASK

MEDICAL

The Task Medical Hi-Lo Couch 3 Sections 1 Motor is a versatile and reliable option for medical settings. Its design allows for easy positioning, with smooth height adjustment controlled by a 360° foot bar system. The couch is designed for durability and patient comfort, featuring a face hole and gas strut-operated sections. Retractable castors and levelling feet ensure stability and mobility.

### Product Features

- Type: 3-Section examination couch
- Upholstery: navy blue
- Mobility: retractable castors and levelling feet
- Patient Comfort: face hole, gas strut operated sections
- Warranty: Frame – 5 years, electrical – 3 years, upholstery – 1 year

CODE  
TK3SECNBFCB



### Product Specifications

- Size (Flat): L1980 x W700 mm
- Adjustable Height: L480 – H900 mm
- Height Adjustment: Single electric actuator with 360° foot bar system
- Back Section Adjustment: 90° to -33°
- Foot Section Adjustment: 41° to -24°
- Weight Capacity: SWL 250kg

Outer Carton



[www.constar.com.au](http://www.constar.com.au)

