

Hi-Lo Exam Couch 3 Sections 3 Motor - Navy Blue

TASK

MEDICAL

The Task Medical Hi-Lo Couch 3 Sections provides precise adjustability with three electric motors controlling the height, headrest, and footrest. Designed for comfort and efficiency, it features a sturdy frame, retractable castors, and levelling feet, making it a reliable choice for various medical environments.

Product Features

- Three electric motors for height, head, and foot adjustments
- Heavy-duty PVC vinyl upholstery – fire, stain, mildew, oil, water, and wear-resistant
- Adjustable back section (90° to -33°) & foot section (41° to -24°)
- Retractable castors & pad levelling feet for stability
- Face hole for enhanced patient comfort

CODE
TK3SEC3MNB



Product Specifications

- Flat size: L 1980 x W 700 mm
- Height adjustability: L 480 – H 900 mm
- Warranty: frame – 5 years, electrical – 3 years, upholstery – 1 year
- Safe Working Load (SWL): 250kg

Outer Carton



www.constar.com.au

