

Hi-Lo Exam Couch with Side Rails 3 Sections 1 Motor - Black

TASK

MEDICAL

The Task Medical Hi-Lo Couch features side rails that swing away for easy access, with upholstered padded vinyl and steel frames. This versatile 3 section bed allows for easy adjustments in height, back, and foot sections. Designed for medical efficiency and patient comfort, the couch ensures reliable service, with a high safe working load and ergonomic electric height adjustments for practitioners.

Product Features

- Swing away side rails for easy patient access
- Heavy-duty PVC vinyl upholstery – fire, stain, mildew, oil, water, and wear-resistant
- Retractable castors & levelling feet for smooth mobility and stability
- Face hole for patient comfort during exams
- Gas strut operated sections for smooth adjustments
- Also available in Navy Blue

CODE
TK3SECIMSRBLK



Product Specifications

- Flat size: L 1980 x W 700 mm
- Height adjustability: L 480 – H 900 mm
- Warranty: frame – 5 years , electrical – 3 years , upholstery – 1 year
- Safe Working Load (SWL): 250kg

Outer Carton



www.constar.com.au

