

## Multipurpose Hi-Lo Treatment Chair - Navy Blue

# TASK

MEDICAL

The Task Medical Multipurpose Hi-Lo Treatment Chair is designed for versatility, providing optimal support for a wide range of medical procedures. Featuring a three-section adjustable design and two powerful motors, it allows seamless height and position adjustments to accommodate both patient comfort and practitioner needs. The integrated armrests offer additional support, while the ergonomic structure ensures stability and ease of use. Ideal for clinics and healthcare facilities, this chair enhances efficiency, adaptability, and patient care in various treatment settings.

### Product Features

- Adjustable backrest (0 to +78°) & foot section (0 to -90°)
- Electrically controlled height & seat positioning
- Heavy-duty PVC vinyl: fire, stain & water-resistant
- Four retractable medical-grade castors for mobility
- Supports up to 250kg

CODE  
TK3SMTCNB



### Product Warranty

- Frame: 5 years
- Electrical components: 3 years
- Upholstery: 1 year

Outer Carton



[www.constar.com.au](http://www.constar.com.au)

