

# Powder Coated Trolley 2 Drawers

# TASK

MEDICAL

The Task Medical Powder Coated Trolley 2 Drawers offers high-grade durability with smooth castors and an outstanding finish. A great alternative to stainless steel models, each trolley comes with a bottom shelf for additional storage.

## Product Features

- Carbon steel powder-coated for durability
- Two soft-closing drawers for extra storage
- Smooth running castors for easy mobility
- 125mm deep shelf for storage
- Top and bottom side rails on all three sides
- Outstanding finish

## Product Dimensions

- (W) 50cm (D) 50cm (H) 90cm

CODE  
TK2DPC



Outer Carton



[www.constar.com.au](http://www.constar.com.au)

