

## SS Trolley Double Four Drawer

**TASK**  
MEDICAL

With four spacious drawers, the Task Medical Double Trolley Four Drawer is designed for healthcare professionals needing extra storage space. Manufactured from hospital-grade 304 stainless steel, it combines strength with ease of movement and access.

### Product Features

- High-grade 304 stainless steel throughout
- Soft-closing drawers
- German-made ball-bearing drawer runners
- 125mm deep drawers
- 95mm smooth-rolling casters for easy mobility
- Locks on front casters
- Top and bottom side rails on all three sides
- Smooth welding finish
- Easy movement

### Product Dimensions

- (L) 80cm (W) 50cm (H) 90cm

CODE  
TK4D80N



Outer Carton



[www.constar.com.au](http://www.constar.com.au)

