

SS Trolley Double No Drawer

TASK
MEDICAL

The Task Medical Stainless Steel Double Trolley offers added storage with two spacious compartments. Made from high-quality 304 stainless steel, it is built for durability and ease of use in busy medical environments.

Product Features

- High grade 304 stainless steel throughout
- 95mm smooth-rolling casters for easy mobility
- Locks on front castors
- Top and bottom side rails on all three sides
- Smooth welding finish
- Easy movement
- Maximum weight limit - 50Kg

Product Dimensions

- (L) 80cm (W) 50cm (H) 90cm

CODE
TKOD80N



Outer Carton



www.constar.com.au

