

SS Trolley No Drawer No Rails

TASK
MEDICAL

The Task Medical Stainless Steel Trolley with No Drawers and No Rails is designed for versatility and durability in medical environments. Constructed from hospital-grade 304 stainless steel, this trolley offers a sleek, easy-to-clean surface and smooth mobility for seamless use in healthcare settings.

Product Features

- High grade 304 stainless steel throughout
- 95mm smooth-rolling casters for easy mobility
- Locks on front castors
- Top and bottom side rails on all three sides
- Smooth welding finish
- Easy movement

Product Dimensions

- (W) 50cm (D) 50cm (H) 90cm

CODE
TKOD50NR



Outer Carton



www.constar.com.au

