

SS Trolley Single Four Drawer with Lock

TASK
MEDICAL

The Task Medical Single Four Drawer Trolley with lock offers secure and organised storage for medical professionals. The top drawer can be locked for added security, while the smooth-rolling casters ensure easy movement within busy medical settings.

Product Features

- Key lock on top drawer for security of contents (spare key included)
- High-grade 304 stainless steel throughout
- Soft-closing magnetic drawers
- German-made ball-bearing drawer runners
- 125mm deep drawers
- 95mm smooth-rolling casters for easy mobility
- Locks on front casters
- Top and bottom side rails on all three sides
- Double safety folds
- Smooth welding finish

Product Dimensions

- (W) 50cm (D) 50cm (H) 90cm

CODE
TK4D50LK



Outer Carton



www.constar.com.au

